

Medical Devices

## DimAir Line

**Intensive Care** 



NIMV



## NIMV HELMET

Single patient use Helmet for Non Invasive Mechanical Ventilation; it can be used also in CPAP and BI-PAP for the weaning from the ventilator, both in adults and children patients. It has 6 ISO connectors placed directly on the rigid ring, to offer more EFFICACY, PRACTICALITY and COMFORT. The version provided with pneumatically sealed zip, allows a quick and easy access to the patient.

**EFFICACY** The connection to the rigid ring reduces the delay of the trigger reaction time. **PRACTICALITY** The double Right and Left symmetrical connectors allow a practical and easy connection with the equipment placed on any side.

COMFORT

The connection between the Helmet and the tubes directly on the rigid ring makes the gas diffusion more comfortable and it reduces the inner noise.

## **TECHNICAL FEATURES**

- Opening zip with pneumatic seal\*
- Single patient use
- Made of biocompatible plastic material
- Soft snd flexible straps with protective pads made of antibacterial CoolMax<sup>®</sup> material to reduce ulcerations
- 6 Right and Left symmetrical connectors
- 2 fixing eyelets on the top of the helmet for an additional/alternative conncection
- 22M, 22F ISO connectors
- Latex Free

**ADVANTAGES** 



The "eyelet" additional connection prevents from ulceration under the ampits and improves inspiratory trigger reaction and the synchrony between patient and ventilator

- Effective, Practical and Comfortable
- Trigger reaction delay < 0,1"</li>
- · Quick acces to the patient's face
- Possibility to do nursing, feeding and weaning operations without removing the helmet
- Gases inlet and outlet can be inverted or connected on the same side; in this way the patients can also be placed on their side

## REFERENCES

)	CODE Standard * ZIP		SIZE	Neck circumference cm.	Internal volume** Lt.
	Standard	ZIP		circumerence cm.	(except head)
	600/5585	600/5586	XXS	23 ÷ 27	9
	600/5228	600/5234	XS	27 ÷ 31	9
	600/5229	600/5235	S	31 ÷ 35	9
	600/5230	600/5236	М	35 ÷ 39	10
	600/5231	600/5237	L	39 ÷ 43	10
	600/5232	600/5238	XL	43 ÷ 47	11
	600/5233	600/5239	XXL	47 ÷ 51	11
	600/5580	600/5582	3XL	≥ 51	11

\*\*+/- dimensional tollerances, excluding the collar and the patient's head

Single pack

Multiple pack with 5 pieces



REG. N. MED 24039 UNI EN ISO 9001:2008 UNI EN ISO 13485:2012 UNI EN ISO 14001:2004



DIMAR

Via G. Galilei, 6 41036 Medolla -MO- I Tel: +39 0535 611336 Fax: +39 0535 611328 Mail: dimar.srl@tiscali.it